


Eye Surgery Stretcher

Maximum Surgical Access.
Simple Functionality.



Stryker's eye surgery stretcher offers the **greatest head-end clearance in the industry**, giving surgeons ample leg room to manoeuvre while performing delicate eye procedures. The firmly-locking work surface and mobile design combine to increase patient security and surgical throughput.



Flexible tube doubles as air delivery and drape support.



Enhanced clearance and dual-articulating headpiece maximises surgical access.



Yellow



Grey



Light Blue



Red



Navy



Purple

Multiple bumper colour options allow you to customise the stretcher for identification purposes.

Enhanced head-end clearance for **supreme surgeon access**

Ideal for either **superior or temporal** approach

Positioning controls tuck neatly within the headpiece and allow full articulation

Optional **contoured or flat head-piece** comfortably supports patient during procedures

Eye Surgery Stretcher



Stryker SA
Cité-Centre
Grand-Rue 90
1820 Montreux
Switzerland

t: +41 21 966 12 01
f: +41 21 966 12 00

www.stryker.eu

Standard Features

- 3" (8 cm) Enhanced Comfort Mattress
- 8" (20 cm) Omni Surface casters with wheel covers
- Dual cylinder pneumatic backrest
- Dual pedestal hydraulics
- Dual side hydraulic controls
- Fold-down siderails
- Four IV receptacles
- Central-locking brakes with foot-end brake & steer pedal
- Hands-free Trendelenburg
- Integrated storage tray with oxygen bottle holder
- Retractable fifth-wheel steering
- Rubber bumpers
- Stationary foot-end litter
- Uni-lower pedal
- Dual-articulation, enhanced clearance headpiece

Specifications

Model Number	1089	CE
Overall Length	90"	(229 cm)
Overall Width (siderails up)	31.5"	(80 cm)
Weight Capacity	500 lbs.	(227 kg)
Height Range (to litter top)		
High	34"	(86.4 cm)
Low	22.25"	(56.5 cm)
Litter Positioning		
Backrest	0° - 90°	
Knee Gatch	0° - 30°	
Trend./Reverse Trend.	±18°	
Headpiece		
Height Adjustment	±4.5"	(11 cm)
Headpiece Tilt	±45°	
Patient Surface	26" x 87"	(66 x 221 cm)
Siderails	13" x 55"	(33 x 140 cm)
Caster Diameter	8"	(20 cm)

Stryker reserves the right to change specifications without notice.

Optional Features

- 4" (10 cm) Ultra Comfort Mattress
- Pre/post-op shoulder pads assembly (crank backrest only)
- Air delivery system/drape support
- Crank backrest
- Crank knee gatch
- Defibrillator tray/foot extender
- Temporal or superior wristrest
- Three-sided hydraulic controls
- Footboard/chartholder
- Three-sided brake & steer pedals
- Integrated patient transfer/arm board
- IV caddy
- Restraint straps
- Serving/instrument tray
- Siderail pads
- Surgery accessory rail
- Three-position siderails
- Upright oxygen bottle holder
- Oxygen bottle retaining ring
- Serving/instrument tray holder
- Head pads - Ultra Comfort
 - 3" (8 cm), flat
 - 3" (8 cm), concave
- IV poles
 - Permanent folding two- or three-stage at foot end
 - Removable

This document is intended solely for the use of healthcare professionals.

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Omni Surface, Stryker. All other trademarks are trademarks of their respective owners or holders.

The products listed above are CE marked according to the Medical Device Directive 93/42/EEC.



Literature Number: MTX6301262005EN Rev F
MTX-RRD/GS 03/12

Copyright © 2012 Stryker