

Cub

Paediatric Crib



Easily-lowered siderails for unobstructed access

Swing-open access doors allow quiet, convenient patient care

Retracting backrest helps eliminate head-end gap

Retractable fifth wheel helps improve mobility

Safety for Smaller Patients

Responsible paediatric care demands a stretcher specifically designed around safety, access and ease of use. The Cub paediatric crib supports safe and efficient care practices while the warm, friendly appearance cultivates a healing environment for small children and their families.

Accent Color Options



Optional protective top equipped with intermediate support position on each side.



Retracting backrest reduces the gap against the siderail while swing-open access doors allow easy patient care.



Lowered siderails provide unobstructed access.



190 Hydraulic Model



190 Fixed Height Model

Standard Features

- 3" (8 cm) vinyl mattress
- Manually-operated retracting backrest
- Manual foot prop
- Siderail with 9" (23 cm) mandatory stop
- Counterbalance system on siderail for easy raising and lowering
- Oxygen bottle holder and retaining ring
- Fixed foot-end rail
- Chart shelf or monitor tray
- Head/Foot end roller bumpers
- IV poles - permanent folding two-stage, head and foot end

Hydraulic Model: All Standard Features Plus

- 6" (15 cm) Omni Surface casters
- Central brake control
- Dual side-mounted foot control hydraulics with uni-lower pedal
- Hands-free Trend./reverse Trend

Optional Features

- 4" (10 cm) Ultra Comfort Mattress
- Siderail without 9" (23 cm) mandatory stop
- Access doors at all four corners
- Fully-lowerable foot-end rail – provides unobstructed access around entire unit when all siderails are down
- Retractable fifth-wheel steering
- IV caddy
- IV poles
 - Permanent folding - head and foot end - three-stage
- Pneumatic assist backrest
- Retractable protective top - availability TBC
- Retractable protective top with intermediate support position - availability TBC
- Upright oxygen bottle holder

Specifications

	190 Hydraulic		190 Fixed Height	
• Model Number				
• Overall Length & Width				
With Accessory Supports	72.25" x 37.5"	(184 x 95 cm)	72.25" x 37.5"	(184 x 95 cm)
• Weight Capacity	400 lbs	(181 kg)	400 lbs	(181 kg)
• Height Range (to litter top)				
High	40"	(102 cm)	33"	(84 cm)
Low	33"	(84 cm)	33"	(84 cm)
• Litter Positioning				
Manual Activated	0° / 25° / 40° / 55°		0° / 25° / 40° / 55°	
Pneumatic Assist	0° - 50°		0° - 50°	
Trend./Reverse Trend.	±9.5°		0° / 6° / 12° with manual foot prop	
Foot Section Angle	0° / 6° / 12°		0° / 6° / 12°	
• Patient Surface	29.5" x 57"	(75 x 145 cm)	29.5" x 57"	(75 x 145 cm)
• Caster Diameter	6"	(15 cm)	5"	(13 cm)

Stryker reserves the right to change specifications without notice.

Warranty

- One year parts, labour and travel

This document is intended solely for the use of healthcare professionals.

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **Cub**, **Omni Surface** and **Stryker**. All other trademarks are trademarks of their respective owners or holders.

The products listed above are CE marked according to the Medical Device Directive 93/42/EEC.



Literature Number: MTX6201262005 Rev. E
MTX6365/GS 03/10

Copyright © 2010 Stryker

Stryker SA
Cité-Centre
Grand-Rue 90
1820 Montreux
Switzerland

t: +41 21 966 12 01
f: +41 21 966 12 00

www.europe.stryker.com