

# Fluoroscopy Stretcher

---

**U-shaped base and full-length C-arm access enhance imaging capabilities**

---

**Fibersin pneumatic backrest at both ends eliminates need to reposition patient**

---

**Siderails fold down and out of the way for improved access**

---

**Optional accessory rail accommodates surgical attachments**

---

**Four-wheel brake system helps secure treatment platform**

---



## **Ease of Imaging**

Access, functionality and safety are the hallmarks of Stryker's fluoroscopy stretcher. This specialty platform features full-length C-arm capability and fold-down siderails for maximum access, while dual-end hydraulic foot controls allow you to safely raise or lower the patient surface.



Dual pneumatic backrest improves imaging access.

# Fluoroscopy Stretcher

stryker®



## Standard Features

- 2" (5 cm) mattress
- 6" (15 cm) non-conductive casters
- Dual end-mounted foot control hydraulics
- Dual fiberesin pneumatic backrest
- Fold-down siderails
- Four corner roller bumpers
- Four IV receptacles
- Four-wheel, steel-ring brake system with dual end activators
- Radiolucent litter top
- Locking-caster steering with dual end activators
- U-shaped base

## Specifications

• <b>Model Number</b>	<b>1080</b>	
• <b>Overall Length</b>	87"	(221 cm)
• <b>Overall Width</b>	30.5"	(77 cm)
• <b>Weight Capacity</b>	500 lbs	(228 kg)
• <b>Height Range (to litter top)</b>		
High	34"	(86 cm)
Low	24"	(61 cm)
• <b>Litter Positioning</b>		
Dual Backrest Operation	0° - 90°	
Trend./Reverse Trend.	±8°	
• <b>Patient Surface</b>	26" x 84"	(66 x 213 cm)
• <b>Siderails</b>	13" x 55"	(33 x 140 cm)
• <b>Fluoroscopy Window</b>	19.5" x 66"	(50 x 168 cm)
• <b>Caster Diameter</b>	6"	(15 cm)

*Stryker reserves the right to change specifications without notice.*

## Optional Features

- Footboard/chartholder
- IV poles
  - Permanent folding two- or three-stage
  - Removable
- Integrated accessory rail
- Restraint straps
- Upright oxygen bottle holder

**This document is intended solely for the use of health-care professionals.**

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **Fluoroscopy, Stryker**. All other trademarks are trademarks of their respective owners or holders.

The products listed above are CE marked according to the Medical Device Directive 93/42/EEC.

Stryker SA  
Cité-Centre  
Grand-Rue 90  
1820 Montreux  
Switzerland

t : +41 21 966 12 01  
f : +41 21 966 12 00

[www.europe.stryker.com](http://www.europe.stryker.com)

Literature Number: MTX6501262005  
MTX5495-09/GS 05/09



Copyright © 2009 Stryker