

# Gynnie

## OB/GYN Stretcher



**Adjustable footrests feature simple, one-handed operation**

**4" thick Ultra Comfort mattress provides unparalleled support**

**Sealed mattress and clean lines prevent fluid collection**

**Compact design with pop-up push handles and retractable fifth wheel improves maneuverability**

### Intuitive Patient Positioning

Gynnie combines the trademark durability and features of a standard Stryker stretcher with supreme OB/GYN functionality. Ideal for pelvic exams and general use, this multi-functional stretcher is ideal for both pelvic exams and general use.



Glide footrests out, position patient and lower foot section to provide excellent pelvic access.



Gynnie can be configured with optional OB/GYN accessories.

## OB/GYN Stretcher



### Standard Features

- 4" (10 cm) Ultra Comfort Mattress with sealed seams
- 8" (20 cm) Omni Surface casters with wheel covers
- 26" (66 cm) wide patient surface
- Dual cylinder pneumatic backrest
- Dual pedestal hydraulics with constant rate descent
- Dual side-mounted foot control hydraulics
- Retractable fifth-wheel steering
- Four-wheel, steel-ring brakes with dual side and head activators
- Glideaway siderails
- Integrated storage tray with oxygen bottle holder
- Permanently-attached folding foot section
- Permanently-attached foot supports
- Pop-up push handles
- Roller bumpers
- Three IV receptacles
- Unique one-handed foot section release lever

### Specifications

• <b>Model Number</b>	<b>1061</b>	
• <b>Overall Length</b>	81.5"	(207 cm)
• <b>Overall Width</b>		
Siderails Up	31.5"	(80 cm)
Siderails Down	26.75"	(68 cm)
• <b>Weight Capacity</b>	500 lbs	(228 kg)
• <b>Height Range (to litter top)</b>		
High	36"	(91 cm)
Low	22.5"	(57 cm)
• <b>Litter Positioning</b>		
Backrest	0° - 90°	
Foot Section	0° - 90°	
Trend./Reverse Trend.	±18°	
• <b>Patient Surface</b>	26" x 77"	(66 x 196 cm)
• <b>Siderails</b>	13" x 39"	(34 x 99 cm)
• <b>Caster Diameter</b>	8"	(20 cm)

Stryker reserves the right to change specifications without notice.

### Warranty

- One year

### Optional Features

- C-spine cassette film holder
- Calf supports and brackets
- Defibrillator tray/footboard/extender
- Backrest X-ray cassette holder
- Hoop for drainage basin
- IV caddy
- IV poles
  - Permanent folding two- or three-stage
  - Tethered
  - Standard removable
- Patient restraints
- Serving/instrument tray
- Siderail pads
- Upright oxygen bottle holder

**This document is intended solely for the use of healthcare professionals.**

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **Glideaway, Gynnie, Omni Surface, Stryker**. All other trademarks are trademarks of their respective owners or holders.

The products listed above are CE marked according to the Medical Device Directive 93/42/EEC.



Literature Number: MTX7001262005  
MTX5494-09/GS 05/09

Copyright © 2009 Stryker

Stryker SA  
Cité-Centre  
Grand-Rue 90  
1820 Montreux  
Switzerland

t : +41 21 966 12 01  
f : +41 21 966 12 00

www.europe.stryker.com

