

Stretcher Chairs



Model 5050

Increase Patient Throughput

Ideal for surgery centers and small spaces, Stryker’s stretcher chairs maximize patient throughput without compromising patient care and outcomes. Choose optimal positioning from horizontal to upright to achieve superior clinical performance and patient comfort.

Stretcher easily articulates from a flat to a chair position

Pneumatic backrest maneuvers 0 to 90 degrees

Dual hydraulic jacks and four-wheel, steel-ring brakes to optimize stability

Foot section can be moved dependently or independently of backrest

Eye stretcher chair model features dual-articulation headpiece for eye surgery



Model 5051

The eye stretcher chair provides all the benefits of the standard stretcher chair plus a dual-articulation, enhanced-clearance headpiece for surgical positioning.

Accent Color Options



Stretcher

Chairs

stryker®



Model 5050



Model 5051



Model 5051 Supine Position

Standard Features

- 2" (5 cm) vinyl mattress
- 6" (15 cm) casters
- 24" (61 cm) wide patient surface
- Adjustable height footrest
- Dual cylinder pneumatic backrest
- Dual pedestal hydraulics
- Fixed position pushbar (model 5050 only)
- Four IV receptacles
- Four-wheel, steel-ring brakes with dual side activators
- Hands-free Trend.
- Swing-down siderails with comfort armrest
- Thermoformed ABS hood with integrated storage tray and oxygen bottle holder

Model 5051: All Standard Features Plus

- Dual-articulation, enhanced-clearance headpiece
- Head pads—Ultra Comfort
 - 3" (8 cm), 2" (5 cm) or 1" (2.5 cm) flat
 - 3" (8 cm) or 2" (5 cm) concave

Optional Features

- 2" (5 cm) or 3" (8 cm) Enhanced Comfort Mattress
- 3" (8 cm) Ultra Comfort Mattress
- Concave head pad (model 5051 only)
- Independent foot section
- IV caddy
- IV poles
 - Tethered
 - Standard removable
- Heel stirrups and brackets
- Leg supports and brackets
- Multi-position pushbar (model 5050 only)
- Integrated patient transfer/arm board
- Restraint straps
- Extended warranty available

Specifications

• Model Numbers	5050/5051	
• Overall Length	76"	(193 cm)
• Width		
Overall	30"	(76 cm)
Litter	27.5"	(70 cm)
• Weight Capacity	400 lbs	(182 kg)
• Height Range (to litter top)		
High	33.5"	(85 cm)
Low	22"	(56 cm)
• Litter Positioning		
Backrest	0° - 90°	
Leg	0° - 80°	
Trend./reverse Trend.	±18°	
• Siderails		
Standard	10" x 31"	(25 x 79 cm)
Transfer	9" x 31"	(23 x 79 cm)
• Patient Surface	24" x 74"	(61 x 188 cm)
• Caster Diameter	6"	(15 cm)

Stryker reserves the right to change specifications without notice.

Model 5051: All Optional Features Plus

- 3" (8 cm) inflatable head pad
- Adjustable arm board and clamp (surgery accessory rail required)
- Air delivery system/drape support
- Surgery accessory rail (not available with patient transfer siderails)
- Temporal or superior wristrest

This document is intended solely for the use of healthcare professionals.

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademark or service mark: Stryker. All other trademarks are trademarks of their respective owners or holders.

The products listed above are CE marked according to the Medical Device Directive 93/42/EEC.



Literature Number: MTXLSS45504C
MTX/GS 03/09

Copyright © 2009 Stryker

Stryker SA
Cité-Centre
Grand-Rue 90
1820 Montreux
Switzerland

t : +41 21 966 12 01
f : +41 21 966 12 00

www.europe.stryker.com