

Stryker Prime Series



The Big Wheel can reduce start-up force up to 50 percent and steering effort up to 60 percent.

Mobile. Safe. Efficient.

Extraordinarily mobile and user-friendly, the Stryker Prime Series stretchers are designed to greatly reduce the physical strain on clinicians while enhancing patient care and comfort. Advanced mobility options such as the Zoom Motorized Drive System or Big Wheel make transporting heavy patients easy, while electric functionality allows patients to adjust their own position without calling a caregiver for assistance. Fully equipped and highly configurable, the Stryker Prime Series allows you to create the ideal solution to meet your specific needs.



Zoom Motorized Drive System on the Stryker Prime stretcher virtually eliminates manual pushing.



Lift Assist Backrest uses a drop-seat design which uses the patient's weight to help raise the head of the bed while cradling the patient to help reduce boosting.



Electric functionality gives optional control to the patient.

Advanced mobility, with optional Big Wheel or Zoom Motorized Drive System, helps make transport safe and efficient

Optional electric-powered base and litter

Pioneer mattress offers superior pressure redistribution and comfort

Power-washable design with 700 lb capacity

Stryker Prime Series

stryker®



1105



1115



1125

Standard Features

- Thermoformed ABS base with extra capacity storage tray
- 700 lb weight capacity
- 26-inch wide patient surface with 3-inch Enhanced Comfort Mattress
- Less than 21-inch low height
- Integrated utility tray
- Litter-mounted oxygen bottle holder
- Power-washable design (1105 and 1115 only)
- Four-wheel, center-locking brake system with dual end activators
- Dual side-mounted foot control hydraulics with uni-lower pedal
- Glideaway Siderails with integrated foot end push/pull handles
- 4 integrated IV receptacles
- +/- Trendelenburg
- 8-inch Omni Surface Casters with integrated wheel covers
- Pneumatic backrest/stationary foot
- Constant descent hydraulics
- Dual end brake/steer controls
- 1105: All Standard Features Plus**
- Fifth Wheel steering
- 1115: All Standard Features Plus**
- Big Wheel steering
- 1125: All Standard Features Plus C-US UL**
- Zoom Motorized Drive System

Specifications

Model Number

26-inch Litter

30-inch Litter

• Model Number	1105 (Prime with Fifth Wheel) 1115 (Prime with Big Wheel) 1125 (Prime with Zoom Drive)	1105 (Prime with Fifth Wheel) 1115 (Prime with Big Wheel) 1125 (Prime with Zoom Drive)
-----------------------	---	---

Overall Length

(1105 and 1115)	85"/216 cm	85"/216 cm
(1125)	86"/218 cm	86"/218 cm

Overall Width

Siderails Up	34"/86.4 cm	38"/96.5 cm
Siderails Down (1105 and 1115)	30.25"/76.8 cm	30.5"/77.5 cm
Siderails Down (1125)	30.25"/76.8 cm	31"/78.7 cm

Weight Capacity

700 lbs/318 kg	700 lbs/318 kg
----------------	----------------

Height Range

High	34"/86.4 cm	34"/86.4 cm
Low (1105 and 1115)	20.75"/53 cm	20.75"/53 cm
Low (1125)	23"/58 cm	23"/58 cm

Litter Positioning

Backrest	0° – 90°	0° – 90°
Knee Gatch	0° – 40°	0° – 40°
Trend./Reverse Trend. (1105 and 1115)	±17°	±17°
Trend./Reverse Trend. (1125)	±16°	±16°

Maximum Angle

87.5°	87.5°
-------	-------

Patient Surface

26" x 75.25"/66 x 191 cm	30" x 75.25"/76 x 191 cm
--------------------------	--------------------------

Siderails

14" x 57.75"/36 x 147 cm	14" x 57.75"/36 x 147 cm
--------------------------	--------------------------

Caster Diameter

8"/20 cm	8"/20 cm
----------	----------

Warranty

- One year parts, labour and travel

Stryker reserves the right to change specifications without notice.

Optional Features

- 30-inch (76 cm) wide patient surface
- 4-inch (10 cm) Enhanced Comfort Mattress
- 4-inch (10 cm) or 5-inch (13 cm) Ultra Comfort Mattress
- Pioneer Pressure Redistribution Mattress
- Four-sided brake and steer controls (1105 and 1115 only)
- Pop-up push handles (1105 and 1115 only)
- Knee Assist No-crank Hydraulic Knee Gatch
- Exclusive Lift Assist Backrest design with drop seat
- Three-sided hydraulic controls (1105 and 1115 only)
- Defibrillator tray

- Defibrillator tray/foot extender/charting surface
- Dual-ended siderail release
- Footboard/chartholder
- IV caddy
- IV poles (permanent, folding two- or three-stage)
- Serving tray
- Serving tray holder/footboard
- Siderail pads
- Upright oxygen bottle holder
- Restraint straps
- Integrated transfer board

Electric Functionality

- Electric option (only on 30 inches wide)
- Electric patient controls
- Foot end nursing controls
- Powered lift (1105 and 1115 only)

This document is intended solely for the use of healthcare professionals.

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Big Wheel, Chaperone, Glideaway, Knee Assist, Lift Assist, Omni Surface, Pioneer, Stryker, Stryker Prime and Zoom. All other trademarks are trademarks of their respective owners or holders.

The products listed above are CE marked according to the Medical Device Directive 93/42/EEC.



Literature Number: MTX369111709
MTX6678/GS 05/10

Copyright © 2010 Stryker



Stryker SA
Cité Centre
Grand-Rue 90
1820 Montreux
Switzerland

t: +41 21 966 12 01
f: +41 21 966 12 00
www.stryker.eu