

Transport Stretcher

Retractable fifth wheel provides improved mobility

Four-sided brake and steer controls offer caregiver flexibility

Fold-down siderails help ensure patient safety

Optional knee gatch provides added patient comfort



Easy Patient Transport

The simple design of Stryker's transport stretcher provides safe and efficient transport for an economical value. This stretcher features fifth wheel mobility for easier traction and cornering, in addition to four-wheel, central locking brakes for optimum stability.



Stryker's transport stretcher is streamlined to provide quick and effective patient transport.

Bumper Color Options



Blue Light Blue Gray Purple Red Yellow



Standard Features

- 3" (8 cm) vinyl mattress on 26" (66 cm) litter
- 3" (8 cm) Enhanced Comfort Mattress on 30" (76 cm) litter
- 8" (20 cm) Omni Surface™ casters
- Fold-down siderails on 26" (66 cm) litter (fold to foot)
- Fold-down siderails on 30" (76 cm) litter (fold to head)
- Four IV receptacles
- Four-wheel, central-locking brakes system with dual-end activators
- Integrated storage tray with oxygen bottle holder
- Retractable fifth-wheel steering
- Roller bumpers
- Stationary foot-end litter
- Dual cylinder pneumatic backrest
- Dual pedestal hydraulics
- Dual side-mounted foot control hydraulics
- Hands-free Trendelenburg positioning
- Wheel covers

Optional Features

- 4" (10 cm) Enhanced Comfort Mattress
- 4" (10 cm) Ultra Comfort Mattress
- 5" (13 cm) Ultra Comfort Redistribution Mattress
- Defibrillator tray
- Defibrillator tray/footboard/extender
- Footboard/chartholder
- Four-sided brake & steer controls
- Pop-up push handles
- IV caddy

Specifications

	26" Litter		30" Litter	
• Model Number	747		748	
• Overall Length	83"	(211 cm)	83"	(211 cm)
• Overall Width	30"	(76 cm)	34"	(86 cm)
• Weight Capacity	500 lbs	(228 kg)	500 lbs	(228 kg)
• Height Range (to litter top)				
High	36"	(91 cm)	36"	(91 cm)
Low	21.5"	(55 cm)	21.5"	(55 cm)
• Litter Positioning				
Backrest	0°- 90°		0°- 90°	
Trend./Reverse Trend.	±18°		±18°	
Knee Gatch			0°- 30°	
• Patient Surface	26" x 75.5"	(66 x 192 cm)	30" x 75.5"	(76 x 192 cm)
• Siderails	13" x 55"	(33 x 140 cm)	13" x 55"	(33 x 140 cm)
• Caster Diameter	8"	(20 cm)	8"	(20 cm)

Stryker reserves the right to change specifications without notice.

Warranty

- One year parts, labour and travel

This document is intended solely for the use of healthcare professionals.

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **Big Wheel, Omni Surface, Stryker**. All other trademarks are trademarks of their respective owners or holders.

The products listed above are CE marked according to the Medical Device Directive 93/42/EEC.



Literature Number: MTX6101252005 Rev. G
MTX6705/GS 06/10

Copyright © 2010 Stryker

Stryker SA
Cité-Centre
Grand-Rue 90
1820 Montreux
Switzerland

t: +41 21 966 12 01
f: +41 21 966 12 00

www.europe.stryker.com