

# Surgistool

The Ideal  
Surgeon's Stool



## The Surgistool customisable seat configurations with fully adjustable backrest and armrests

provide optimal support to meet your specific positioning preference. Superior pedal access and a stable yet mobile design combine to deliver the best possible clinical performance.



Teardrop Seat



Wedge Seat



Saddle Seat  
(optional)

**Three-spindle base** provides superior control and pedal access

**Dual-density foam seat and adjustable backrest** offer maximum support

**Optional padded armrests** allow versatile positioning and stability

**Multiple seat configurations** offer personalised comfort

# Surgistool



**Teardrop Seat**



**Wedge Seat**



**Saddle Seat (optional)**

## Standard Features

- Adjustable contoured backrest
- Dual-density foam seat cushions (except saddle seat option)
- Foot-controlled, hydraulic, three-spindle base
- Front end locking system
- Thermoformed plastic hood
- Teardrop or wedge seat

## Optional Features

- 2-inch (5 cm) riser kit
- Saddle seat
- Multiposition armrest

## Specifications

Model Number	830	
Weight Capacity	300 lb.	(136 kg)
Height Range		
High	28.5"	(72 cm)
Low	21"	(53 cm)
Height Range with Riser Kit		
High	30.5"	(77 cm)
Low	23"	(58 cm)
Seat Cushion Style	Teardrop	
	Wedge	
	Saddle	

*Stryker reserves the right to change specifications without notice.*

### This document is intended solely for the use of healthcare professionals.

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademark or service mark: **Surgistool, Stryker**. All other trademarks are trademarks of their respective owners or holders.

The products listed above are CE marked according to the Medical Device Directive 93/42/EEC.



Literature Number: MTX6601262005EN Rev D  
MTX-RRD/GS 03/12

Copyright © 2012 Stryker